

Aprotic, Polar and Hydrophobic Solvent
For Chemical Reactions and Extractions

MTHP

Chemical Name : 4-Methyltetrahydropyran

Synonym : 4-MeTHP

CAS No. : 4717-96-8

Kuraray

Isoprene chemicals Div.

2507-1

kuraray



Applications of MTHP

- Solvent for Chemical Reactions
- Solvent for Extractions
- Solvent for Polymerizations
- Solvent for Coatings

Chemical Reactions in MTHP

Organometallic Reactions

| | Reference # |
|---|-------------|
| Grignard | 1), 4) |
| Lithiation of aryl bromide | 6) |
| Anionic ortho-Fries rearrangement of carbamates with Lithium compound | 6) |
| Alkylation of aldehyde with dialkyl zinc | 6) |
| LAH Reduction | 1) |

Transition Metal Catalyzed Reactions

| | |
|---|--------|
| Suzuki – Miyaura Coupling | 7), 9) |
| Sonogashira-Hagiwara cross coupling | 1), 7) |
| Palladium-Catalyzed three component coupling | 8) |
| Nickel-Catalyzed Direct Arylation of Aliphatic Amides | 2) |
| Platinum-Catalyzed cyclodimerization of alkynes | 8) |
| Ruthenium-Catalyzed Hydrogenation of carbonyl compounds | 5) |
| Olefin Metathesis | 3), 7) |

Substitute for halogenated solvent

| | |
|--------------------------------|--------|
| TEMPO catalyzed Oxidation | 1), 4) |
| Dess-Martin Oxidation | 1) |
| Swern Oxidation | 1) |
| Epoxidation of allylic alcohol | 1) |
| Amidation of Acid | 1) |
| Ring-Closing-Metathesis | 1) |

Others

| | |
|---|----------|
| Tin or Silicon radical-mediated addition and reduction | 1) |
| Michael addition | 1) |
| Wittig | 1) |
| Dehydration condensation - Esterification and Acetalization | 4) |
| Reduction | |
| Halogenation, Dehalogenation | |
| Protection & Deprotection for Peptide synthesis | 10), 11) |

See Slide 17 for the references

Features of MTHP & Benefits with MTHP in Chemical reactions

Features

- ✓ Aprotic & Polar
- ✓ Hydrophobic

- ✓ B.P. is higher than THF
- ✓ Stable against Acids etc.
- ✓ Better Reaction yield depending on reaction


Benefits

- Good in balance of polarity/solvency and extraction yield
- Quick and clear phase separation with water
- No need to change the solvent throughout the entire process, reaction and washing/extraction
- Easy to recycle MTHP
- Less COD/BOD in wastewater
- Reduced CO₂ emission
- Higher in Reaction Rate
- Less by-products in Reaction & Quench.
→ Higher purity of target products
- Lower Raw Material Cost


Physical Properties

Solvents are arranged in order of increase of “Solubility of solvent into water”.

| | B.P. [°C] | M.P. [°C] | Density [20°C] | Viscosity [mPa·s] | F.P. [°C] | Solubility of solvent into Water [wt%] | Solubility of Water into Solvent [wt%] | Azeotrope with Water [°C], Water wt% | Log Pow |
|---------|--------------|--------------|-------------------|----------------------|--------------|---|---|--|------------|
| Toluene | 111 | -95 | 0.87 | 0.56 | 4.4 | 0.05 | 0.03 | 85, 20 | 2.73 |
| CPME | 106 | -140 | 0.86 | 0.55 | -1 | 1.1 | 0.3 | 83, 16 | 1.59 |
| DCM | 40 | -97 | 1.32 | 0.44 | --- | 1.3 | 0.2 | 40, 5 | 1.25 |
| MTHP | 105 | -92 | 0.86 | 0.78 | 6.5 | 1.5 | 1.4 | 85, 19 | 1.90 |
| EA | 77 | -84 | 0.90 | 0.43 | -4 | 8.7 | 3.3 | 70, 8 | 0.73 |
| 2MTHF | 80 | -136 | 0.85 | 0.6 | -11 | 14 | 4.4 | 71, 11 | 0.77 |
| THF | 65 | -109 | 0.89 | 0.55 | -15 | ∞ | ∞ | 64, 6 | 0.46 |
| Water | 100 | 0 | 1 | | --- | ∞ | ∞ | --- | --- |



Polarity and Compatibility with water



Despite its high polarity, MTHP also exhibits high hydrophobicity.

The unique property of MTHP are presumed to be derived from its six-membered cyclic ether structure.

Source : Shoji Kobayashi et al., Chem. Asian J. **2019**, 14, 3921 – 3937

Phase separation with water -1



<Test Method >

Added 15ml of water colored green to each test tube.

Added 15ml of each of different solvents to each test tube.

The test tubes were set to the holder, shaken together for about 20 sec. and left to stand for about 3 hours at 23°C / 73°F.

THF : No phase separation with water.

2MTHF and Ethyl Acetate:
The level of water layer was elevated.

MTHP : No change in the level of water layer.

Phase separation with water -2

Length of time after the shaking
for the emulsion layer to disappear

| Solvent | Time [Seconds] | | |
|---------------|--------------------------------|---------|-----------------------------|
| | 23°C | 40°C | 40°C |
| | No Salt | No Salt | With NaCl 20wt% to Water |
| THF | No Phase Separation with Water | | |
| 2MTHF | 30 | 20 | 10 |
| MTHP | 5 | 2 | <1 |
| CPME | 2 | | |
| Toluene | 8 | | |
| Ethyl Acetate | 16 | | |
| DCM | >60 | | |


<Test Method >

Added 15ml of water colored green to each test tube.

Added 15ml of each of different solvents to each test tube.


The test tubes were set to the holder, shaken together for about 20 sec.
and left to stand at 23°C / 73°F.

The time was measured for the emulsion layer to disappear and
For the interface level stabilize.



Click

Solubility of Organic & Inorganic salts into each solvent at 22°C



Compatibility of the solvents
with Organic salts & Inorganic salts


THF >> 2MTHF \doteq MTHP > CPME

Stability of solvents

Autoxidation


< Test condition >

- Under air at 25°C, without stabilizer.
- POV was measured by Iodometry method.



Residual rate of NBL(n-butyl lithium)


50°C



Acidic condition

< Test condition >


- 20%HCl/Solvent=1:1(wt.) @50°C




Basic condition

< Test condition >

- 20%NaOH/Solvent=1/1(wt.) @50°C




Chemical Reaction : Grignard



| Solvent | Yield [%] | | |
|---------|---------------|----------------|-------|
| | In Org. layer | In Aqua. Layer | Total |
| THF | 0.0 | 0.0 | 0.0 |
| 2MTHF | 70.5 | 2.5 | 73.0 |
| MTHP | 83.9 | 0.9 | 84.9 |

MeCl(b.p. -24°C), MeBr(b.p. 4°C) are not easy to handle because of their low b.p. though MeI(b.p. 42°C) is high enough in b.p. to handle.



However, in the reaction of methyl iodide and magnesium in THF, the methyl Grignard reagent cannot be prepared in a high yield due to side reactions.



The side reaction was inhibited, and target compound was produced in good yield with MTHP.

Chemical Reaction : Lithiation and Boronization


- Preparation of Raw material for Suzuki-Miyaura coupling -



Yield of boron compound were good no matter which solvent is used when the reaction was conducted at -78°C , Which is too low to conduct on a commercial scale.

Only MTHP system had good yield of boronic acid
When the reaction was conducted at -20°C , Which is feasible temperature on a commercial scale.

Chemical Reaction : Suzuki – Miyaura Coupling




| Solvent | Reaction Temp. * [°C] | Reaction time [hrs] | Yield Of coupling compounds |
|-------------|--------------------------|------------------------|-----------------------------------|
| THF | 62 | 8 | 70 |
| 2MTHF | 70 | 4 | 82 |
| MTHP | 82 | 2 | 93 |
| CPME | 82 | 2 | 89 |
| Toluene | 84 | 2 | 88 |

MTHP system was good in both reaction rate and Yield as Toluene system was, which is a popular solvent for this reaction.

* Azeotropic reflux with Water, a co-solvent

Chemical Reaction : TEMPO catalyzed Oxidation

- Substitute for DCM -



| Solvent | Yield [%] | | | Notes |
|---------|-----------------|----------|-------------------------|--------------------|
| | Carboxylic Acid | Aldehyde | Alcohol Raw Material | |
| DCM | 75 | 0.9 | 0.6 | Low in Selectivity |
| CPME | 64 | 15 | 13 | Low in conversion |
| MTHP | 89 | 8.5 | 0.9 | |

MTHP can be a substitute for DCM, Dichloromethane, which is a hazardous solvent.

General Information

- Product Name : MTHP
- CAS No. : 4717-96-8
- Package : Drum (Net 170kg)

Specifications










- Appearance : Colorless transparent
- Purity (%GC) : ≥ 99.0
- Water (ppm) : ≤ 200

* Contains Stabilizer : BHT 20ppm

Regulatory status

| Government Inventory list | Status | Remarks |
|------------------------------|-----------------------------|---|
| ENCS (Japan) | Present | |
| EU-REACH (Europe) | Registered | Registrant : Kuraray Europe GmbH (Importer) Tonnage Band 10-100 MT |
| TSCA (USA) | Registered | SNUR |
| IECSC (China) | Provisionally Registered | Tonnage Band : 1-10 MT |
| K-REACH (Korea) | | OR registration Tonnage Band : 10-100 MT |

Hazardous Information

| | MTHP | THF | 2MTHF |
|-----------------------------------|--|---|---|
| Pictogram |    |    |    |
| Flammable liquids | Category 2 | Category 2 | Category 2 |
| Skin corrosion/irritation | Category 1 | | Category 2 |
| Serious eye damage/eye irritation | Category 1 | Category 2 | Category 1 |
| Skin sensitization | Category 1 | | |
| Acute toxicity (oral) | | Category 4 | Category 4 |
| Carcinogenicity | | Category 2 | |
| Specific target organ toxicity | | Category 3 | |
| Reference | SDS, Kuraray Dec 13, 2023 | SDS, Thermo Fisher Scientific Dec 6, 2024 | SDS, Thermo Fisher Scientific Sep 22, 2023 |

Dissolution of Resins in the different solvents

| RESIN | MTHP | THF | 2MTHF | CPME |
|---------------|------|-----|-------|------|
| ABS | D | D | D | D |
| Acrylic resin | D | D | D | D |
| Polystyrene | D | D | D | D |
| PVC | D | D | D | P |
| Polycarbonate | P | D | D | P |
| Fluoro rubber | P | P | D | P |
| SBR | P | P | D | P |
| Polyurethane | P | P | D | P |
| Butyl rubber | P | ND | P | ND |
| Nylon-6 | ND | ND | ND | ND |
| Phenol resin | ND | ND | ND | ND |
| PP | ND | ND | ND | ND |
| HDPE | ND | ND | ND | ND |
| LDPE | ND | ND | ND | ND |
| Teflon | ND | ND | ND | ND |

< Test Method >

Test panels of the resin were dipped into each solvent at 50°C for 7days. After the dipping, weight of the test panels were measured and compared with original one.

D : Dissolved completely

P : Dissolved 30% to less than 100%

ND : Dissolved 0% - 30%

Reference

- 1) Shoji Kobayashi, Araki Masuyama et al., *Chemistry An Asian Journal*, **2019**, Volume 14, Issue 21, 3921-3937
- 2) Miki Iyanaga, Yoshinori Aihara, and Naoto Chatani, *The Journal of Organic Chemistry*, **2014**, 79, 11933-11939
- 3) Anna Kajetanowicz, Karol Grela et al., *ACS Sustainable Chem. Eng.*, **2020**, 8, 49, 18215-18223
- 4) Yuusuke Saitou, *Speciality Chemicals Magazine*, **2023**, May/June, 30-32
- 5) Yoichi Hoshimoto et al., *Tetrahedron Chem*, **2024**, 9, 100059
- 6) Marco Blangetti et al., *Eur. J. Org. Chem.* **2024**, 27, e202400313
- 7) Bhaskar N. Thorat et al., *Organic Process Research & Development*, **2022**, 26, 3, 480-492
- 8) Mamoru Tobisu et al., *Organometallics*, **2020**, 39, 1678-1682
- 9) Haruki Mizoguchi, Akira Sakakura et al., *Chem. Sci.*, **2022**, 13, 9580-9585
- 10) Daniel Sejer Pedersen et al., *RSC Adv.*, **2020**, 10, 42457-42492
- 11) Hisata, Y., Washio, T., Takizawa, S. et al., *Nat Commun.*, **2024**, 15, 3708

Contact Information

Kuraray America, Inc.

Advances Chemicals
3700 Bay Area Blvd. Suite
#680, Houston, TX 77058

U.S.A.

P +1 281 283 1790

chemical.sales@kuraray.com

Kuraray Co., Ltd.

Isoprene Chemicals division
Tokiwabashi Tower
2-6-4, Otemachi, Chiyoda-ku
Tokyo 100-0004

Japan

<https://www.isoprene-chemicals.com/index.html>

Kuraray Trading (Shanghai) Co., Ltd.

Unit 2106, 2 Grand Gateway, 3 Hongqiao
Road, Xuhui District, Shanghai 200030,

China

P +86 21 6407 9182

<https://www.kuraray-trading.co.jp/eng/>

Kuraray Europe GmbH

BU Chemicals
Philipp-Reis-Straße 4
65795 Hattersheim am Main

Germany

P +49 69 305 35844

chemicals@kuraray.com

Kuraray India Private Ltd.

Prius Platinum 2nd Floor B Wing,
D3 District
Center Saket New Delhi 110017,

India

P +91 11 4090 4400

<https://kurarayindia.co.in/>

Kuraray Trading (Taipei) Co., Ltd.

4 F., No. 372, Linsen N. Rd.,
Zhongshan Dist.,
Taipei City 104

Taiwan (R.O.C.)

P +886 2-2521-8831

<https://www.kuraray-trading.co.jp/eng/>

© Kuraray Co., Ltd. 2023

Brand/Product name, Brand/Product name and Brand/Product name is a/are trademark(s) or registered trademark(s) of Kuraray Co., Ltd. or its affiliates. Trademarks may not be applied for or registered in all countries. The information, recommendations and details given in this document have been compiled with care and to our best knowledge and belief. They do not entail an assurance of properties above and beyond the product specification. The user of our products is responsible for ensuring that the product is suitable for the intended use and conforms to all relevant regulations. Kuraray Co., Ltd. and its affiliates do not accept any guarantee or liability for any errors, inaccuracies or omissions in this document.

kuraray